

PC Algorithm Application

Neill Stanford, CHP

www.stanforddosimetry.com





Goals/features

- Stand-alone PC based dose algorithm
- Calculate doses using standard algorithm
- Allow adjustments to element responses
- Perform “what-if” testing using simulated element responses
- Modular coding to allow efficient update of response matrix or algorithm design
- Formal reporting to printer and export to spreadsheet application.




General design

- Two modes of operation
 - Responses – input element responses (keyboard or file)
 - Doses – input doses and responses are estimated
- Matrix calculation
 - Specify fields to use for fit
 - Specify number of iterations
- Three reporting options including export to file
- On-line algorithm documentation

Options

- Background subtraction
- Position-specific correction factors
- Units
- Matrix calculation



Responses mode

- Input element readings, with background and correction factors if desired.
- The dose is calculated using the standard algorithm.

InLight NVLAP 12/15/03 - [Enter TLD Response Data]

File Setup Operation Help

New Group... Calculate Details Print Report... Exit

Enter a comment for this TLD response analysis here.

Net Element Responses

E1	E2	E3	E4
790	380	360	260

Background Values

0	0	0	0
---	---	---	---

Correction Factor

1	1	1	1
---	---	---	---

Multiplier Divisor

Algorithm Processing Results :

Photon (keV)	82.2	Shallow Dose	Eye Dose	Deep Dose	mrem
Reported Dose		1176	196	203	beta
Matrix dose		1215	231	234	2000 TESTS

Show Matrix

Doses mode –

Tests performance of algorithm under ideal conditions without the cost, time and added bias of irradiations.

1. Select a field
2. Specify a dose level
3. Program calculates and sums responses based on inputs and response matrix
4. Doses are calculated as normal, with the percent difference displayed for each dose.

The screenshot shows the 'InLight NVLAP 12/15/03 - [Enter Dose Data]' window. It features a menu bar (File, Setup, Operation, Help) and a toolbar with buttons for 'New Group...', 'Calculate', 'Details', 'Print', 'Report...', and 'Exit'. A list of radiation sources is on the left, with 'Cs-137 662keV' selected. A 'Shallow Dose in mrem' input field contains '100'. The 'Total Shallow mrem' field shows '1200.0'. Below this is a 'Summary of Entered Doses (mrem)' table:

TI-204	1000
Am-241 60keV	100
Cs-137 662keV	100

At the bottom, a 'Calculated Element Responses (mR*)' table is displayed:

789.48	379.62	360.37	260.72
--------	--------	--------	--------

Algorithm Processing Results:

Photon (keV)	82.4	Shallow Dose	Eye Dose	Deep Dose	mrem
Reported Dose		1175	196	204	beta
Delivered Dose		1200	208	210	
% Difference		-2.1%	-5.7%	-2.9%	
Matrix dose		1210	236	236	200 TESTS

4

Any mixture can be input, in this case a large beta component with both x-ray and Cs photons added.

Full page detailed report

Detailed report provides all of the interim calculations to allow troubleshooting results.


1. Corrected element responses
2. Photon dose estimates
3. Non-photon calculations
4. Final doses

OSL Results					
<i>InLight NVLAP 12/15/03</i>					
<i>Full Page Format</i>					
Description	Testing alg performance			Printed	May 13, 2004
Run Time	05-13-2004 15:06:04			Run By	CNP
		Units	mrem		
Background	0.00	0.00	0.00	0.00	Path option: Normal
Correction factors	1.00	1.00	1.00	1.00	
<hr/>					
1					
Corrected element responses					
1	<u>E1 (OW)</u>	<u>E2(PL)</u>	<u>E3(A1)</u>	<u>E4(Cu)</u>	<u>R34 (E3/E4)</u>
	789.50	379.60	360.40	260.70	1.38
Initial calculations					
2	<u>Photon dose calculations</u>	<i>Based on E2</i>	<u>Shallow</u>	<u>Deep</u>	<u>Photon energy (keV)</u>
		<i>Based on E4</i>	186.40	206.52	82.38
			178.72	201.09	
3	<u>Non-photon dose calculations</u>				
	Non photon on E1	408.78	<u>NetE1/NetE2</u>	100.00	
	Non photon on E2	-1.74			
	<u>Calculated neutron</u>		0.00	<u>Calculated shallow beta</u>	988.34
	<i>Neutron correction factor =</i>			<u>Calculated lens of eye beta</u>	0.00
4	Reported Doses				
	<u>mrem</u>	<u>Shallow</u>	<u>Eye</u>	<u>Deep</u>	<i>Flags</i>
	Beta	988.34	0.00		
	Neutron	0.00	0.00	0.00	
	Photon	186.40	196.00	203.81	82.4 keV
	Totals	1,174.75	196.00	203.81	mrem beta




File input/output

- Spreadsheet files can be used for input, and results can be exported back to spreadsheet format.



On-line documentation

Full algorithm documentation with response tables, correction factors and flow charts is available through the help menu.





STANFORD DOSIMETRY

Dosimetry program support since 1988

- Algorithm development
- Program documentation
- Technical support
- Training workshops

Email: info@stanforddosimetry.com

website www.stanforddosimetry.com

other www.svmargarita.com

Mail: 619 Commercial Ave. #28
Anacortes, WA 98221

